


MANUAL Dishwashing Method

Méthode de lavage MANUEL de la vaisselle



Instructions

Step/Étape 3

Sanitize

for at least 45 seconds, in:

1. hot water, not lower than 77°C (170°F), OR

2. water, not lower than 24°C (75°F) with:

• iodine solution, 25 mg/L (25 ppm), OR

• chlorine solution, 100 mg/L (100 ppm), OR

• quaternary ammonium solution, 200 mg/L (200 ppm).

Note: Test kit for sanitizer must be available.

Désinfecter

pendant au moins 45 secondes :

1. eau chaude, minimum de 77 °C (170 °F), OU

2. eau, minimum de 24 °C (75 °F) contenant :

• une solution d'iode, 25 mg/l (25 ppm), OU

• une solution de chlore, 100 mg/l (100 ppm), OU

• une solution d'ammonium quaternaire, 200 mg/l (200 ppm).

Veuillez noter : une trousse d'analyse pour la concentration de désinfectant doit être disponible.